

Don't Let Smoke Get in Your Eyes!

The Geiger Smoke Evacuator is designed to be the most cost effective system available on the market.

Its rugged design and quality construction ensures years of operation. Geiger model SA-1A smoke evacuator includes unit and one pack of tubing. Optional wand, footswitch and hose stand are available. Filter packs are available in three different configurations: FIL3 (for non-smoke procedures such as wrinkle removal, podiatry and dental procedures). 2-year warranty.



Specifications Include:

- Particulate Filtration: 99.999% at 0.2 micron (worst case)
- Odor absorption via activated charcoal for approx. 60 minutes of use
- Filter Longevity: approx. 30-45 minutes of smoke suction
- Airflow Control: approx. 0.5 CFM to Max
- Noise Level: 60-65 db at 1 meter
- Weight: 18 pounds
- Voltage: 110V or 220V available
- Manufactured per FDA and UL guidelines

Order by contacting your Geiger Medical Technologies Distributor or Geiger Medical Technologies.

ACCESSORIES



GEIGER SMOKE EVACUATOR

Dimensions: 11" W x 12" H x 10" D
Item # SA-1A

FILTER PACK WITH TUBING

Item # FIL1, FIL2 or FIL3

OBEDIENT HOSE STAND™

Item # SA-H1

9" ALUMINUM WAND

Item # WND1

www.geigermedical.com